

Catalogue

一、	Overview	1
	1. Purposes	
	Study the validity period of kits and samples under different storage conditions.	
	2. Product name and specification	1
	3. Reference product composition and product performance indicators	
	3.1 Reference product composition	
	3.2 Performance	
_,	Stability study	2
	2.1Shipping Stability (Unopened kit)	
	2.2. Real-time stability study (Unopened kit)	
	2.3. In-use/Opened Kit Stability - Sample diluent	
	2.4. In-use/Opened Kit Stability - Test card	

一、Overview

1. Purposes

Study the validity period of kits and samples under different storage conditions.

2. Product name and specification

Novel coronavirus (2019-nCoV) antigen detection kit (colloidal gold method), referred to as "this kit". Specification: 20test/box

3. Reference product composition and product performance indicators

3.1 Reference product composition

Table 1 Reference product composition

Types	No.	Microorganism category
	QT01	Staphylococcus aureus
	QT02	Streptococcus pneumoniae
	QT03	Measles virus
	QT04	Mumps virus
	QT05	Adenovirus type 3
	QT06	Mycoplasma pneumoniae
	QT07	Parainfluenza type 2
	QT08	Metapneumovirus
	QT09	Coronavirus OC43
Negative reference	QT10	Coronavirus 229E
riegative reference	QT11	Bacillus parapertussis
	QT12	Influenza B virus (Victoria line)
	QT13	Influenza B virus (Yamagata series)
	QT14	H1N1 (2009) influenza virus
	QT15	Influenza A H3N2 virus
	QT16	Epstein-Barr virus
	QT17	Enterovirus CA16
	QT18	Rhinovirus
	QY01	SARS-CoV-2
	QY02	SARS-CoV-2
Positive reference	QY03	SARS-CoV-2
	QY04	SARS-CoV-2
	QY05	SARS-CoV-2
D. C.C.	CV01	SARS-CoV-2 (Moderately Positive)
Repetitive reference	CV02	SARS-CoV-2 (Weakly positive)
Lowest detection limit reference	S01	SARS-CoV-2

3.2 Performance

Samples to be tested:

Number of samples (Number of repeats): 18 negative enterprise reference samples (external negative controls), 5 positive enterprise reference samples (external positive controls), 2 repeatability enterprise reference samples (external positive controls) and 4 LoD enterprise reference samples (external positive controls). A known positive clinical sample diluted to 3x LoD by negative nasopharyngeal swab specimen, a known positive clinical sample diluted to 5x LoD by negative nasopharyngeal swab specimen, a known negative clinical sample.

Source: from clinically diagnosis SARS-CoV-2 cases and non SARS-CoV-2 cases.

Storage: store at -70° C.

Number of replicates: 1 replicate per negative enterprise reference sample, positive enterprise reference sample and LoD enterprise reference sample. 10 replicates per repeatability enterprise reference sample. 5 replicates per 3x LoD clinical sample, 5x LoD clinical sample and negative clinical sample.

Acceptance criteria:

3.2.1 Appearance

Use visual inspection with normal vision or corrected vision under natural light: the components of the kit should be complete, and the outer packaging should be clean and free of contamination. The test card is sealed and stored without damage and pollution. The sample diluent is colorless, clear and homogeneous, without visible impurities or precipitates.

3.2.2 The width of membrane strip

Take one film strip randomly from the kit and measure with vernier caliper. The width of film strip shall not be less than 2.5mm.

3.2.3 Liquid migration speed

Select three test cards randomly and place them horizontally. Measure the distance L (mm) between the center of the sampling hole and the distal end of the observation window. Use the sample diluent as the sample to be tested and record the time t (s) required for the liquid to move to the far end of the observation window with a stopwatch. Calculate L / T (mm / min), which is the migration speed. The migration speed of liquid should not be less than 10 mm / min.

3.2.**4** Coincidence rate of positive reference

The positive enterprise reference samples should be tested once, and the results should be positive.

3.2.5 Coincidence rate of negative reference

The negative enterprise reference samples should be tested once, and the results should be negative.

3.2.6 Limit of detection (LoD)

Dilute the LoD enterprise reference sample into S1, S2, S3 and S4 by 2 times, 4 times, 8 times and 16 times, respectively. Test each sample once, S1 and S2 results should be positive, and S3 and S4 results can be negative or positive.

3.2.7 Repeatability

The results of 10 tests on each repeatability enterprise reference sample should be positive and uniform in color.

3.2.8 Accuracy of clinical sample test

Test two clinical positive samples and one clinical negative sample, and the test results should be consistent with the known results.

二、Stability study

2.1Shipping Stability (Unopened kit)

• Claims: $40 \, \text{C}$ for less than $30 \, \text{days}$.

• *Number of test points:* testing at 0,10,20,30,35 days.

• Baseline: 0 day from bottling

• *Product Lots to be tested:* 20200801, 20200802, 20200803

• Experimental protocol: Test the kits before transportation. Then the test kits are placed in a shaker that can maintain a constant temperature of 40 $^{\circ}$ C for transportation simulation, which is shaken at a certain frequency to simulate transportation conditions. The kits are tested at 0,10,20,30,35 days.

• **Results:** Test results are as follows.

Lot number	Performand Tim		Day 0	Day 10	Day 20	Day 30	Day 35
	Appear	rance	Qualified	Qualified	Qualified	Qualified	Qualified
	The width membran (mn	n of the ne strip	4.03	4.10	4.04	4.06	4.08
	7 · · · · · · ·	,•	26.7,	26.45,	26.32,	25.89,	26.32,
	Liquid mi	_	26.71,	26.74,	26.51,	25.78,	26.28,
	speed (mi	m/min)	26.35	26.62	26.49	26.11	26.42
	Coincidence positive re	•	5/5 (+)	5/5 (+)	5/5 (+)	5/5 (+)	5/5 (+)
20200001	Coincidence rate of negative reference		18/18 (-)	18/18 (-)	18/18 (-)	18/18 (-)	18/18 (-)
20200801	Repeatab	CV01	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)
	ility	CV02	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)
		S1	+	+	+	+	+
	1 D	S2	+	+	+	+	+
	LoD	<i>S3</i>	+	+	+	+	+
		<i>S4</i>	-	+	+	ı	+
	Accuracy	3*LoD	+	+	+	+	+
	of	5*LoD	+	+	+	+	+
	clinical sample test	Negati ve	-	-	-	-	-
	Appear	ance	Qualified	Qualified	Qualified	Qualified	Qualified
	The width membran (mn	ne strip	4.05	4.07	4.03	4.02	4.06
	7 · · · · · · ·	,•	25.74,	25.86,	26.12,	26.31,	26.45,
	Liquid mi	_	25.94,	25.59,	26.11,	26.33,	26.51,
20200802	speed (mi	m/min)	26.7	25.65	26.21	26.41	26.84
	Coincidence positive re	_	5/5 (+)	5/5 (+)	5/5 (+)	5/5 (+)	5/5 (+)
	Coincidence negative re	-	18/18 (-)	18/18 (-)	18/18 (-)	18/18 (-)	18/18 (-)
	Repeatab	CV01	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)
	ility	CV02	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)

		C 1					
		S1	+	+	+	+	+
	LoD	S2	+	+	+	+	+
		<i>S3</i>	+	+	+	+	+
		S4	+	+	+	-	+
	Accuracy	3*LoD	+	+	+	+	+
	of	5*LoD	+	+	+	+	+
	clinical sample test	Negati ve	-	-	-	-	-
	Appear	ance	Qualified	Qualified	Qualified	Qualified	Qualified
	The width membran (mn	e strip	3.92	3.89	3.94	3.97	3.91
	Liquid mi speed (mi	_	26.05, 25.83, 25.94	26.21, 26.45, 25.98	26.06, 25.86, 25.81	26.11, 26.2, 26.13	25.97, 25.92, 26.02
	Coincidence rate of positive reference		5/5 (+)	5/5 (+)	5/5 (+)	5/5 (+)	5/5 (+)
20200803	Coincidence negative re		18/18 (-)	18/18 (-)	18/18 (-)	18/18 (-)	18/18 (-)
20200803	Repeatab	CV01	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)
	ility	CV02	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)
		S1	+	+	+	+	+
	<i>I</i> D	S2	+	+	+	+	+
	LoD	<i>S3</i>	-	+	+	+	+
		S4	-	-	+	+	-
	Accuracy	3*LoD	+	+	+	+	+
	of	5*LoD	+	+	+	+	+
	clinical sample test	Negati ve	-	-	-	-	-

• Conclusions: The results show that the test kit stored for 35 days under simulated transportation environment performs normally. Under the principle of ensuring the validity of the test kit, we recommend that the transportation condition of the test kit be stored at 40 $^{\circ}$ C, less than 30 days.

2.2. Real-time stability study (Unopened kit)

•Claims: Store at 2 $\mathcal{C}\sim 30$ \mathcal{C} , valid for 18 months.

•Number of test points: testing at 0,1,2,3,6,12,18,20 months

• Baseline: 0 day from bottling

•Product Lots to be tested: 20200801, 20200802, 20200803

•Experimental protocol: Test the kits as soon as they are bottled. After the kits were qualified, each batch of kits were divided into two parts, one was stored at 2 ~ 8 °C and the other was stored at 30 °C for real-time stability study. The kits were stored for 1,2,3,6,12,18,20 months and then taken out for testing. Real time stability is still in progress.

• Results: Test results are as follows.

Baseline Result

Performance Lot no	umber	20200801	20200802	20200803
Аррес	arance	Qualified	Qualified	Qualified
	membrane strip m)	4.04	4.06	4.07
Liquid migration	n speed (mm/min)	26.24, 26.51, 25.77	25.98, 25.94, 25.88	26.01, 26.11, 26.08
	rate of positive rence	5/5 (+)	5/5 (+)	5/5 (+)
	ate of negative rence	18/18 (-)	18/18 (-)	18/18 (-)
Domontaliita	CV01	10/10 (+)	10/10 (+)	10/10 (+)
Repeatability	CV02	10/10 (+)	10/10 (+)	10/10 (+)
	S1	+	+	+
LoD	S2	+	+	+
LoD	S3	+	+	+
	S4	-	+	-
Accuracy of	3*LoD	+	+	+
clinical sample	5*LoD	+	+	+
test	Negative	-	-	-

Results of kits stored at 2~8 $^{\circ}$ C

Lot numb er	Performan	Time	1 month	2 months	3 months	6 months	12 months	18 months	20 months
	Appeara	ınce	Qualified	Qualified	Qualified	On-going			
	The widt the memb strip (m	rane	3.91	4.05	4.01				
	Liquie	d	26.19,	25.74,	26.53,				
	migration		26.13,	26.05,	25.6,				
	speed (mm/min)		25.94	24.95	26.13				
20200 801	Coincidence rate of positive reference		5/5 (+)	5/5 (+)	5/5 (+)				
	Coincide rate of neg referen	gative	18/18 (-)	18/18 (-)	18/18 (-)				
	Repea (CV01	10/10 (+)	10/10 (+)	10/10 (+)				
	tabilit y	CV02	10/10 (+)	10/10 (+)	10/10 (+)				
	LoD	S1	+	+	+	-		-	

		<i>S</i> 2	+	+	+			
		<i>S3</i>	+	-	-			
		S4	-	-	+			
	Accur	3*Lo						
	acy of	D	+	+	+			
	clinic	5*Lo						
	al	D	+	+	+			
	sampl	Negati						
	e test	ve	-					
	Appea	ırance	Qualified	Qualified	Qualified	On-going		
	The w	idth of						
	the men	mbrane	4.02	3.99	4.01			
		(mm)						
		ruid	26.13,	26.35,	26.05,			
	migr		26.35,	25.94,	24.95,			
		nm/min)	25.83	24.95	25.94			
		idence 		- (- ()	- (- ()			
		positive	5/5 (+)	5/5 (+)	5/5 (+)			
		rence						
		idence	10/10 ()	10/10 ()	10/10 ()			
		negative	18/18 (-)	18/18 (-)	18/18 (-)			
20200 802		ence CV01	10/10 (+)	10/10 (+)	10/10 (+)			
802	Repea tabilit	CVOI	10/10 (+)	10/10 (+)	10/10 (+)			
		CV02	10/10 (+)	10/10 (+)	10/10 (+)			
	У	S1	+	+	+			
		S2	+	+	+			
	LoD	S3	-	+	-			
		S4	_	-	+			
	Accur	3*Lo						
	acy of	D	+	+	+			
	clinic	5*Lo						
	al	D	+	+	+			
	sampl	Negati						
	e test	ve	-	-	-			
	Аррес	ırance	Qualified	Qualified	Qualified	On-going		
	The w	idth of				-		
	the me	mbrane	4.03	4.02	4.01			
	strip	(mm)					 	
	Liq	uid	26.05,	26.05,	26.13,		 	
20200	migr	ation	25.94,	25.83,	24.95,			
803	speed (n	nm/min)	26.21	24.95	25.94			
	Coinc	idence						
		positive	5/5 (+)	5/5 (+)	5/5 (+)			
		rence						
		idence	18/18 (-)	18/18 (-)	18/18 (-)			
	rate of 1	negative						

refer	rence					
Repea	CV01	10/10 (+)	10/10 (+)	10/10 (+)		
tabilit y	CV02	10/10 (+)	10/10 (+)	10/10 (+)		
	S1	+	+	+		
LoD	S2	+	+	+		
LoD	<i>S3</i>	•	ı	+		
	<i>S4</i>	+	ı	ı		
Accur acy of	3*Lo D	+	+	+		
clinic al	5*Lo D	+	+	+		
sampl e test	Negati ve	-	-	-		

Results of kits stored at 30 $^{\circ}\mathrm{C}$

Lot numb er	Perform	Time ance	1 month	2 months	3 months	6 months	12 months	18 months	20 months
	Appe	earance	Qualifi ed	Qualified	Qualified	On-going			
	The width of the membrane strip (mm)		4.05	4.04	4.03				
	_	migration (mm/min)	25.74, 26.21, 26.35	25.42, 25.74, 24.83	26.0, 24.83, 26.53				
	of p	dence rate ositive erence	5/5 (+)	5/5 (+)	5/5 (+)				
2020	of n	dence rate egative erence	18/18	18/18 (-)	18/18 (-)				
0801	Repeat	CV01	10/10 (+)	10/10 (+)	10/10 (+)				
	ability	CV02	10/10 (+)	10/10 (+)	10/10 (+)				
		S1	+	+	+				
		S2	+	+	+				
	LoD	<i>S3</i>	+	-	+				
		S4	-	+	-				
	Accur	3*LoD	+	+	+				
	acy of	5*LoD	+	+	+				
	clinica l sampl e test		-	-	-				
2020	Appe	earance	Qualifi	Qualified	Qualified	On-going			

0802			ed					
0002	The wi	idth of the						
		rane strip	3.97	3.94	3.95			
		mm)	5.57	5.71	5.75			
			26.71,	25.94,	26.35,			
	_	migration	26.05,	26.19,	25.74,			
	speed	(mm/min)	25.6	24.95	25.27			
	Coinci	dence rate						
	of p	ositive	5/5 (+)	5/5 (+)	5/5 (+)			
		erence						
		dence rate						
	of n	egative	18/18	18/18 (-)	18/18 (-)			
		erence	(-)					
			10/10	70/70 ()	70/70 ()			
	Repeat	CV01	(+)	10/10 (+)	10/10 (+)			
	ability	CV02	10/10	10/10 (+)	10/10 (+)			
		C V 02	(+)	10/10 (+)	10/10 (+)			
		S1	+	+	+			
	LoD	S2	+	+	+			
	LoD	<i>S3</i>	-	+				
		S4	-	-	+			
	Accur	3*LoD	+	+	+			
	acy of	5*LoD	+	+	+			
	clinica							
	l	l npl Negative						
	sampl		-	-	-			
	e test							
	4		Qualifi	0 1:0: 1	0 1:0: 1	0 :		
	Appe	earance	ed	Qualified	Qualified	On-going		
	The wi	idth of the						
	memb	rane strip	4.04	4.02	4.03			
	(1	mm)						
	Liquid	migration	25.74,	24.95,	25.6,			
	_	(mm/min)	25.83,	25.94,	26.0,			
			25.42	26.05	24.95			
		dence rate						
		ositive	5/5 (+)	5/5 (+)	5/5 (+)			
2020		erence						
0803		dence rate	18/18					
		egative	(-)	18/18 (-)	18/18 (-)			
	refe	erence						
	_	CV01	10/10	10/10 (+)	10/10 (+)			
	Repeat	5,01	(+)	10,10 (1)	10,10 (1)			
	ability	CV02	10/10	10/10 (+)	10/10 (+)			
			(+)					
		S1	+	+	+			
	LoD	S2	+	+	+			
		<i>S3</i>	-	-	+			
		S4	+	-	-			

Accur	3*LoD	+	+	+		
acy of	5*LoD	+	+	+		
clinica						
l	M					
sampl	Negative	-	-	-		
e test						

2.3. In-use/Opened Kit Stability - Sample diluent

•Claims: Store at 2 $C\sim30$ C, valid for 18 months.

•Number of test points: testing at 1,2,3,6,12,18,20 months

• Baseline: Opening the bottle for the first time, 0 day from bottling.

•Product Lots to be tested: 20200801

•Experimental protocol: Take out a number of test kits and use them for the first time (all bottles of sample diluent are opened to test the physical properties of the kit, the coincidence rate of positive reference, the coincidence rate of negative reference, repeatability, limit of detection and accuracy of clinical sample test). After opening the bottle, the sample diluent is immediately capped and divided into two parts: Part A is kept at 2 ~ 8 °C and part B is kept at 30 °C. The remaining test cards were stored at room temperature (30 °C). In the bottle opening test after 1, 2, 3, 6, 12, 18 and 20 months of storage, 10 drops of sample diluent were vertically extruded and discarded according to the operation in the manual. The appearance of the sample treatment solution was observed (whether it was colorless, clear and homogeneous liquid, without visible impurities or precipitates). Besides the appearance, the positive reference coincidence rate, negative reference coincidence rate, repeatability, limit of detection and accuracy of clinical sample test were also tested. The test is still in progress.

• Results: Test results are as follows.

Baseline Result

4	The width	Liquid	Coincide nce rate			Repeatabilit y		LoD				Accuracy of clinical sample test		
Appear ance	of the membrane strip (mm)	migration speed (mm/min)	of positive reference	of negative reference	CV0	CV0 2	S 1	S 2	S 3	S 4	3*Lo D	5*Lo D	Negati ve	
Qualifi ed	4.05	26.53, 26.13, 25.94	5/5 (+)	18/18 (-)	10/1 0 (+)	10/1 0 (+)	+	+	-	-	+	+	-	

Appearance Results of Sample diluent

Temp $^{\mathcal{C}}$	2~8 °C								
Time	1 month	2 months	3 months	6 months	12 months	18 months	20 months		
Appearance	Qualified	Qualified	Qualified	On-going					
Temp ${\mathcal C}$				30 °C					
Appearance	1 month	2 months	3 months	6 months	12 months	18 months	20 months		
Appearance	Qualified	Qualified	Qualified	On-going					

Results of Coincidence rate of positive reference, Coincidence rate of negative reference,

Repeatability, LoD, Accuracy of clinical sample test

				2~8 °C				
Performance Time		1 month	2 months	3 months	6 months	12 months	18 months	20 months
Coincidence rate of positive reference		5/5 (+)	5/5 (+)	5/5 (+)				
Coincidence rate of negative reference		18/18 (-)	18/18 (-)	18/18 (-)				
Repeatabil	CV01	10/10 (+)	10/10 (+)	10/10 (+)				
ity	CV02	10/10 (+)	10/10 (+)	10/10 (+)				
	S1	+	+	+				
LoD	S2	+	+	+				
	<i>S3</i>	-	+	-				
	S4	+	-	-				
Accuracy	3*LoD	+	+	+				
of clinical	5*LoD	+	+	+				
sample test	Negative	-	-	-				
				30 °C				
Performance Time		1 month	2 months	3 months	6 months	12 months	18 months	20 months
	Coincidence rate of positive reference		5/5 (+)	5/5 (+)				
Coincidence rate of negative reference		18/18 (-)	18/18 (-)	18/18 (-)				
Repeatabil	CV01	10/10 (+)	10/10 (+)	10/10 (+)				
ity	CV02	10/10 (+)	10/10 (+)	10/10 (+)				
	S1	+	+	+				
LoD	S2	+	+	+				
	S3	+	-	-				
	S4	-	-	+				
Accuracy	3*LoD	+	+	+				
of clinical	5*LoD	+	+	+				
sample test	Negative	-	-	-				

2.4. In-use/Opened Kit Stability - Test card

- Claims: Store at 30 C, stable for 1 hour
- *Number of test points:* testing at 0,20,40,60,80 minutes
- Baseline: Tearing off the foil pouch for the first time, 0 day from foil pouch sealing
- Product Lots to be tested: 20200801

- Experimental protocol: Take out a number of test kits, take out the test card from the kit, and use it for the first time (tear the aluminum foil pouch, test the physical properties of the kit, the coincidence rate of positive reference and negative reference, repeatability, limit of detection, and accuracy of clinical sample test), and place the unsealed test card at room temperature (30°C), and store the sample diluent at 2 ~ 8°C. After opening, the test cards were placed for 20 minutes, 40 minutes, 60 minutes and 80 minutes respectively for physical properties, positive reference coincidence rate, negative reference coincidence rate, limit of detection, repeatability and accuracy of clinical sample test.
- Results: Test results are as follows.

· Test results are as jouons.								
Time		0 minute	20	40	60	80		
Performance		0 minute	minutes	minutes	minutes	minutes		
Appearance		Qualified	Qualified	Qualified	Qualified	Qualified		
The width of the membrane strip (mm)		4.06	4.05	4.03	4.04	4.03		
Liquid migration speed (mm/min)		25.68,	25.74,	26.71,	26.35,	26.04,		
		25.76,	26.71,	26.35,	26.21,	26.08,		
		25.94	26.0	27.29	26.53	26.11		
Coincidence rate of positive reference		5/5 (+)	5/5 (+)	5/5 (+)	5/5 (+)	5/5 (+)		
Coincidence rate of negative reference		18/18 (-)	18/18 (-)	18/18 (-)	18/18 (-)	18/18 (-)		
Repeatability	CV01	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)		
	CV02	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)	10/10 (+)		
LoD	S1	+	+	+	+	+		
	S2	+	+	+	+	+		
	<i>S3</i>	+	-	-	+	+		
	S4	+	+	-	-	+		
Accuracy of	3*LoD	+	+	+	+	+		
clinical	5*LoD	+	+	+	+	+		
sample test Negative		-	-	-	-	-		

• Conclusions: The results showed that the performance of test card remained stable for 80 minutes after foil pouch being teared off. Under the principle of ensuring the validity of the test kit, we recommend that test cards should be used as soon as possible within 1 hour after opening the foil pouch.